



Virudet[®] N

A neutral low-foam cleaner for washer disinfector

Technical Information

Description & Application:

Virudet N is for use in washer/disinfectors, AER (Automated Endoscope Reprocessor), other automated decontamination systems and manual cleaning. It is a neutral, concentrated detergent with a wide material compatibility which gently removes blood, proteins, mucus and tissue residues.

This chemistry has been specially developed for the gentle cleaning of delicate surgical instruments, utensils, anaesthesia and endoscope material, allowing simultaneous reprocessing of instruments made of different materials. It has optimal cleaning performance at temperatures of 55-80°C but for endoscopes temperature should not exceed 55°C and can be used for the following: stainless steel instruments (scalpels, tweezers etc.), articulated instruments (scissors, clamps etc.), ward utensils (bowls, feeding bottles etc.), anaesthesia utensils (breathing balloons, tubes etc.), shoes of the surgical team, transport containers made of stainless steel as well as aluminium and endoscopes.

Primary features attributable to Virudet N are:

- Neutral low-foam detergent cleaner;
- Suitable for removal of contaminants from sensitive instrument materials ;
- Uses a blend of complexing and sequestering agents, detergents and corrosion inhibitors that are effective in hard and soft water areas;
- Reduced instrument damage and repair:
- Excellent material compatibility;
- Non-fixative;
- Prevents re-deposition of contaminants and water mark staining;
- Non-corrosive, will not damage delicate coatings.
- Cost effective – improved staff productivity, more flexible;
- Can be used with tap water (saves DI/RO water cost)
- Can be used without neutralizer
- Good cleaning power even at low dosage



**Infection
Controlled**

Applications

Virudet N is ideal as a cleaner in washer/disinfector equipment.

Virudet N is suitable for automated reprocessing of surgical instruments and anaesthesia material, including:

- Stainless steel and articulated instruments (scalpels, tweezers, scissors, clamps, etc.)
- Ward utensils (bowls, feeding bottles etc.)
- Anaesthesia utensils (breathing balloons, breathing tubes etc.)
- Operating Theatre shoes of the surgical team
- Transport containers made of stainless steel as well as aluminium
- Flexible and rigid endoscopes

Compatibility

Virudet N is a neutral cleaner with a broad material compatibility.

However, if there are any concerns then parts should be tested for compatibility prior to processing.

Health and Safety

Virudet N is biodegradable.

The product is classed as non-hazardous in diluted forms, however contact with skin and eyes should be avoided.

See Material Safety Data Sheets for details.

Typical Properties

Appearance	Pale brown liquid
Odour	Odourless
pH,	7 units (neat solution)
Density, g/cm ³ @ 20 °C	1.03 (neat)
Solubility in water	Complete
Flash Point (Abel closed cup)	None

Shelf Life: Concentrate – Not less than 3 years in unopened original containers when stored between 5 °C and 30 °C out of direct sunlight.

Directions for Use

Virudet N is supplied as a concentrated product and should be diluted prior to use.

For Manual & Ultrasonic Cleaning of instruments and equipment:

For cleaning applications, dilute 1 part Virudet N in 100 parts water (i.e. a 1% solution; 10ml/L).

Soak parts for between 5 and 15 minutes in the diluted Virudet N solution at between 50°C to 60°C (allow longer times up to 15 minutes when dealing with dried on proteinaceous deposits).

Note: soak times may vary according to the temperature used and the nature of the contaminants.

After soaking, rinse off thoroughly with clean water (temperature > 40°C) or transfer to the next cleaning stage e.g. manual cleaning, washer/disinfector or ultrasonic equipment.

For Automated Cleaning of instruments and equipment:

Injectors should be adjusted to deliver at 5 - 10ml per litre of water dependant on the expected soil loading of the contaminated equipment and the temperature of the water.

Phase	Water type	Time	Temp	Solution	dosage
Prewash	Cold tap water	3 min	-	-	-
Wash	Hot tap water	10 min	55°C	Virudet N	5 - 10ml/L (0.5 -1%)
Rinse	Hot water	1 min	-	-	-
Rinse	DI water	1 min	-	-	-
Thermo disinfection	DI water	3 min	93°C	Virudet NDA	0.5ml/L or 0.05%
Drying		15 min	110°C		

For thermal disinfection please follow the actual countries rules

Drying time and temp according to the machine manufacturer

General Notes:

Virudet N is not a disinfectant and parts should proceed to disinfection or sterilisation by appropriate methods immediately following cleaning.

Virudet N should not be used at a temperature less than 55°C in washer/disinfector machines.

Neutralisation following cleaning is not required if the above post cleaning rinse cycle is employed.

For maximum effectiveness in Ultrasonic or other soaking applications, the Virudet N solution should be replaced at least daily or when the observed cleaning action falls off.

Do not add non recommended cleaning agents to the Virudet N solution as they can seriously impair performance of the cleaner.

For increase the cleaning power, Virudet N can be safely used in combination with enzymatic boosters in automated reprocessing equipment with improved performance over many products. For details refer to Duo System data sheet.

Dosage over 35°C prevents foam problems.

CONTACT DETAILS:

For UK and Rest of the World:

Amity International,
Libra House, West Street
Worsbrough Dale,
BARNSELEY
S YORKS, S70 5PG
ENGLAND

Tel: +44 (0) 1226 770787

Fax: +44 (0) 1226 770757

For North America:

Amity International,
PO Box 5254,
1704 Denver Road,
ANDERSON,
SOUTH CAROLINA,
SC29623, USA.

Tel: 864 622 2233

Fax: 864 622 2234

E-mail: sales@amityinternational.com

Web site: <http://www.amityinternational.com>

For any further information, please contact your distributor or Amity.

In the event of any technical queries, please contact:
Mr. Ram Singh at the UK/ROW address or by e-mail to:

rsingh@amityinternational.com

Packaging Details:

Virudet N is available in the following packaging.

- 1 L container (available on request)
- 5 L container
- 10 L container
- 25 L container (available on request)
- 200 L drum (available on request)
- 1,000 L IBC (available on request)